



CHAETOMIUM MEDIUM (DSMZ-Medium 188)

NaNO ₃	2.00	g
MgSO ₄ x 7 H ₂ O	0.50	g
KCl	0.50	g
Fe ₂ (SO ₄) ₃ x H ₂ O	0.01	g
KH ₂ PO ₄	0.14	g
K ₂ HPO ₄	1.20	g
Agar	15.00	g
Yeast extract	0.02	g
Distilled water	1000.0	ml

Adjust pH to 7.2. Place a strip of sterile filter paper on cooled agar and inoculate on filter strip.